

POWERFUL 2.0

get back, faster



FAST RECOVERY

COLD & HOT THERAPY UNIT



NO ICE / NO WATER

Therapeutic cooling technology without ice or water inside.

CLOSED & CLEAN SYSTEM

ZAMAR SYSTEM preloaded with glycol doesn't need to refill or drain off liquid.

INNOVATIVE COOLING

Advance technology based on real-time temperature parameters.

STAY FOCUSED ON TEMPERATURE

Set exact temperature you want to reach.

INSPIRED BY YOUR SHAPE

Complete range of Anatomic wraps designed for horse.

Lorenzo De Luca





POWERFUL & STRONGER



powerful than ever

PLUG 'N' PLAY RECOVERY

No set up time and preparation needed.

FRIENDLY USE

All treatment parameters manageable in one screen.

REPEAT YOUR FAVORITE

Possibility to store your personal treatments for frequent use



POWERFUL 2.0



ZAMAR MEDICAL d.o.o.

Sv. Martin, 6 - 52450 - Vrsar - (HR)

Tel: +385 (0)52 496 111 - VAT: 96642432252

info@zamar.care

Features

APPLICATION	COLD - HOT
Temperature Range	-5°C +40°C (32°F - 104°F)
SCREEN	LCD 5" color Touch Screen
WRAPS	Liquid
Glycol	NON TOX freeze
Setup Time Limit	200 h.
Ambient operating	+10°C +30°C (50°F to 95°F)
Salvable therapies	12
Dimensions	cm. 58 x 28 x 30 inc. 23 x 11 x 12
WEIGHT	27 kg. (59 lb.)

- LCD 5 inc, TFT color touch screen
- Easy software, Intuitive Menu icons
- Possibility to select temperature and treatment time
- Set programs that doctors can customize and save
- Totally automated cycle management
- Fast connection/detachment of wraps and pipes with connectors
- Anatomic Wraps pre charged with glycol
- NON-TOX cooling liquid

ONLY FOR VETERINARY USE

Under the FEI regulation and
Veterinary regulations article 1053:
The cooling temperature below 0
degrees C is not permitted.

3ETA081

www.zamar.care

ISO 9001
BUREAU VERITAS
Certification

