

**Z-ONE<sup>®</sup>**

**get back, faster**



# FAST RECOVERY

## COLD & COMPRESSION THERAPY UNIT



### **NO ICE / NO WATER**

Therapeutic cooling technology without ice or water inside.

### **CLOSED & CLEAN SYSTEM**

ZAMAR DRIPLESS SYSTEM preloaded with glycol doesn't need to refill or drain off liquid.

### **INNOVATIVE COOLING**

Advance technology based on real-time temperature parameters.

### **STAY FOCUSED ON TEMPERATURE**

Set exact temperature you want to reach.

### **INTERMITTENT COMPRESSION**

Improves benefits of cold therapy and healing with 3 levels of air-compression.

### **INSPIRED BY YOUR SHAPE**

Complete range of Anatomic wraps designed on your body.

# LIGHTER & STRONGER



**Z-ONE**





# LIGHTER & STRONGER

## Z-ONE



**small and powerful than ever**

### PLUG 'N' PLAY RECOVERY

No set up time and preparation needed.

### FRIENDLY USE

All treatment parameters manageable in one screen.

### ALWAYS ON THE MOVE

Light and handy can easily be with you anytime, anywhere.

### REPEAT YOUR FAVORITE

Possibility to store your personal treatments for frequent use



# Z-ONE®



## Features

<b>APPLICATION</b>	COLD
<b>TEMPERATURE RANGE</b>	+5°C +15°C (41°F 60°F)
<b>SCREEN</b>	LCD 5" color TouchScreen
<b>WRAPS</b>	Liquid / Air Wraps
<b>Glycol</b>	NON TOX freeze
<b>Air Compression levels</b>	mmHg: 50 - 75 - 110
<b>setup time limit</b>	10 h.
<b>Medical class for Wraps</b>	BF
<b>Salvable therapies</b>	1 + 2
<b>DIMENSIONS</b>	cm. 29 x 26 x 30 inc. 11 x 10 x 11
<b>WEIGHT</b>	9 kg. - (20 lb.)

- LCD 5 inc, TFT color touch screen
- Easy software, Intuitive Menu icons
- Possibility to select temperature and treatment time
- Set programs that doctors can customize and save
- Totally automated cycle management
- Fast connection/detachment of wraps and pipes with Dripless connectors
- Anatomic Wraps pre charge with glycol
- NON-TOX cooling liquid

The ZAMAR logo consists of the word "ZAMAR" in a bold, black, sans-serif font, enclosed within a blue oval shape.

**ZAMAR MEDICAL d.o.o.**

Sv. Martin, 6 - 52450 - Vrsar - (HR)  
Tel: +385 (0)52 496 111 - VAT: 96642432252  
info@zamar.care

[www.zamar.care](http://www.zamar.care)

CE 1370

ISO 13485  
BUREAU VERITAS  
Certification



ISO 9001  
BUREAU VERITAS  
Certification

